

Neater Eater Robot



Overview

The **Neater Eater Robotic** enables you to feed yourself. It is simple to use; controlled by a Tablet, plug-in switches or Voice Control.

Tablet menus and help guides demonstrate many options such as setting your mouth position, automatic spoon-wiping and plate-turning. There are fork and snack-holder modes and a routine to scrape the plate at the end of your meal.

- A** On/Off switch
- B** Removable touchscreen Tablet with onscreen help and videos
- C** Robot arm has wide range of movement to reach your mouth
- D** Choice of plug-on cutlery
- E** Choice of plates that rotate

PTO for easy setup!


Please Scan the QR code for the full Instructions.



Quick start guide

- 1 Turn ROBOT on
- 2 Turn TABLET on and follow onscreen instructions
- 😊 Enjoy your meal

(Remember to charge the Robot and Tablet every night)



Removable touchscreen


Tablet for easy access. Replace here between meals to recharge its battery wirelessly.

1 Robot power button.
Press on —powers on and **lights up**, hold to switch off.

2 Tablet power and screen wake/sleep button.

Red Light = charging
(flashing = battery low)
Green Light = fully charged

Plug power cable in here overnight to charge batteries (can be plugged in when eating).



Orange Light = Robot battery charging
Red Light = Power to Tablet charger
Green Light = External power plugged in
Blue Light = Tablet Bluetooth connected

Problem?

Ensure robot and tablet are charged. Hold Tablet power and the volume up buttons until it vibrates. Hold Robot power button until light goes out. Press Robot power button again—it lights up. Wait for Neater App to start.

Plug-in optional switches e.g. to control with foot or head.

Press this question mark icon on screen for help and hold for videos.

The Tablet will guide you through set-up and use.

PLEASE READ ALL SAFETY INSTRUCTIONS.

For further information, video hints, and tips go to www.neater.co.uk

